


Lupus Anticoagulant Profile Testing Algorithm



*Weak lupus anticoagulant may not inhibit DRVVT Mix

** Additional assays may be performed if further clarification or confirmation is necessary based on preliminary review. These may include:

- Coagulation Factor AssaysCoagulation Factor Inhibitor Screen
- Bethesda Units
 Staclot Lupus Anticoagulant
- FibrinogenD-Dimer
- Soluble Fibrin Monomer
 von Willebrand Factor Assays

*** Unfractionated/low-molecular weight heparin or direct thrombin inhibitor (eg, dabigatran, argatroban) warfarin or low concentration of direct anti-Xa inhibitor (eg, rivaroxaban, apixaban, or edoxban)

An interpretive report is included with this profile. If any tests are abnormal, all results will be reviewed by a coagulation consultant.

