

Targeted Genes and Methodology Details for Porphyria Comprehensive Gene Panel

Next-generation sequencing (NGS) and/or Sanger sequencing is performed to test for the presence of variants in coding regions and intron/exon boundaries of the genes analyzed. NGS and/or a polymerase chain reaction (PCR)-based quantitative method is performed to test for the presence of deletions and duplications in the genes analyzed. Confirmation of select reportable variants may be performed by alternate methodologies based on internal laboratory criteria.

Genomic Build: GRCh37 (hg19) unless otherwise specified

Gene	Reference Transcript ^a
ALAD	NM_000031.6
ALAS2	NM_000032.5
CLPX	NM_006660.5
CPOX	NM_000097.7
FECH	NM_000140.4
GATA1	NM_002049.4
HFE	NM_000410.3
HMBS	NM_000190.4
PPOX	NM_000309.5
UROD	NM_000374.5
UROS	NM_000375.3

^a Reference transcript numbers may be updated due to transcript re-versioning. Always refer to the final patient report for gene transcript information referenced at the time of testing.

To verify if a specific region/exon/variant is covered by this assay, contact a laboratory genetic counselor at 800-533-1710.